COUNTER HEIGHT CHAIR WITH CUSHION SEAT

ASSEMBLY INSTRUCTIONS



8PCS 8PCS

2PCS

18PCS

18PCS

2PCS

12PCS

ASSEMBLY TIPS:

- 1. Remove hardware from box and sort by size
- 2. Please check to see that all hardware and parts are present prior to start of assembly.
- 3. Please follow attached instruction in the same sequence as numbered to assure fast & easy assembly.

Warning!

1. Don't attempt to repair or modify parts that are broken or defective. Please contact the store immediately.

2. This product is for home use only and not intended for commercial establishment

PARTS IDENTIFICATION

ARTS IDENTIFICATION				HARDWARE IDENTIFICATION			
A	HEALT	CHAIR BACK	2PCS	1	01	LONG BOLT	8
B		CUSHION SEAT WITH FRAME	2PCS	2	01	MEDIUM BOLT	8
				3	0	SHORT BOLT	2
С	(LEFT FRONT LEG	2PCS	4	0	LOCK WASHER	1
D		RIGHT FRONT LEG	2PCS	5	\bigcirc	FLAT WASHER	1
E		RIGHT SIDE STRETCHER	2PCS	6		ALLEN WENCH	2
F		LEFT SIDE STRETCHER	2PCS	7	0	SCREW	1

FRONT STRETCHER 2PCS

NOTE : QUANTITIES SHOWN ARE FOR TWO CHAIRS STEP 2

STEP 1

G

1.1 Align & hand tighten cushion seat with frame to chair back


bottom rail using short bolt, lock washer & flat washer.

1.2 Align & hand tighten cushion seat with frame to chair back leg

using long bolt, lock washer & flat washer.

1.3 Tighten all bolts with Allen wrench.

1.4 Insert left and right side stretchers to chair back legs.



STEP 3

Complete.



2.2 Then align & fasten left and right front legs to cushion seat with frame using medium bolt, lock washer & flat washer.

TIGNTEN WITH ALLEN WRENCH.

2.3 Fasten front, left & right side stretchers to back, left & right front legs using screw.TIGNTEN WITH SCREW DRIVER (NOT INCLUDED).

2.1 Insert front, left & right side stretchers to left and right front legs.

